

SwallowSTRONG WiFi Guide

The SwallowSTRONG uses a wireless internet connection (WiFi) to sync data from therapy to the SwallowSTRONG Management System (SSMS). If you have access to a WiFi connection, setting it up with your SwallowSTRONG device is quick and easy. Follow the steps below to establish the connection:



Figure 1. Home Screen. Access the WiFi Center by tapping the icon (see above). The Back and Home buttons appear in the same position on every screen.

Step 1: Power the device on and wait until the home screen appears. See Figure 1.

Step 2: Tap the **'WiFi Center'** icon. See Figure 1.

a: A list of available networks appears.

Step 3: Tap the name of the network you would like to join.

a: If the network is password protected, a prompt asking for the password will appear.

b: Enter the password, and tap connect.

c: NOTE: If you don't know the password, contact the network administrator.

Step 4: Once connected, the network will have a **green checkmark** in the column marked **'Status.'**

Step 5: Tap the **back arrow or home screen button** at the bottom of the page to return to the home screen. See Figure 1. From here tap the **'Swallow' icon** to start the SwallowSTRONG program.